




**At the heart of life.  
At the heart of precision.**



# MyVal-1 Study 30 days Outcome



**LOW**<sup>\*</sup>

Device-Related Mortality  
New Permanent Pacemaker  
Incidence of Stroke

**97% Device Success**<sup>\*</sup>

# Myval THV: Designed for Precision in Outcomes


## Hybrid Honeycomb Cell Design Concept

Open cells on upper half to ensure un-jailing of coronary ostia



Closed cells on lower half for high radial strength

≈ 53%  
Frame height 17-21 mm  
≈ 47%



Cobalt alloy frame for high radial strength & radiopacity

Bovine pericardium tri-leaflet valve

External PET buffering to minimize para-valvular leaks

Internal PET sealing cuff for lower profile & puncture resistance

Ø – 20mm, 23mm, 26mm, 29mm<sup>#</sup> - Standard sizes  
Ø – 21.5mm, 24.5mm, 27.5mm - Intermediate Sizes  
Ø – 30.5 mm, 32 mm\* - XL Sizes





Myval THV has been indigenously developed by Meril Life Sciences Pvt. Ltd. & is CDSCO approved.

# - Sizes available now





\* - Contact Meril Representative for information and availability in your country.

Myval THV diameter sizes 30.5 mm and 32 mm are currently not CE approved. Myval THV diameter size 30.5 mm is currently pending approval from CDSCO.

# Myval THV: Size Matrix

| Myval THV<br>Size Matrix<br>& Technical<br>Specifications | Area 314 mm <sup>2</sup><br>17.50 mm<br><br>20 mm | Area 415 mm <sup>2</sup><br>18 mm<br><br>23 mm | Area 531 mm <sup>2</sup><br>19 mm<br><br>26 mm | Area 661 mm <sup>2</sup><br>20.50 mm<br><br>29 mm |
|---|--|---|---|--|
| Perimeter   | 62.83 mm   | 72.26 mm  | 81.68 mm  | 91.11 mm   |
| Native annulus area                                       | 270 - 330 mm <sup>2</sup>  | 360 - 440 mm <sup>2</sup>   | 460 - 560 mm <sup>2</sup>   | 570 - 700 mm <sup>2</sup>  |
| Area-derived diameter                                     | 18.5 - 20.5 mm   | 21.4 - 23.7 mm  | 24.2 - 26.7 mm  | 26.9 - 29.9 mm   |
| Native annulus size by TEE                                | 16 - 19 mm   | 18 - 22 mm  | 21 - 25 mm  | 24 - 28 mm   |

# Myval THV: Additional\* Size Matrix

| Myval THV Size Matrix & Technical Specifications | Myval Intermediate Size  |   |  | Myval XL Size  |  |
|--|--|---|--|--|--|
|  | Area 363 mm <sup>2</sup><br>18.50 mm<br><br>21.5 mm | Area 471 mm <sup>2</sup><br>18.90 mm<br><br>24.5 mm | Area 594 mm <sup>2</sup><br>19.40 mm<br><br>27.5 mm | Area 731 mm <sup>2</sup><br>21.05 mm<br><br>30.5 mm | Area 804 mm <sup>2</sup><br>21.30 mm<br><br>32 mm |
| Perimeter  | 67.54 mm   | 76.97 mm  | 86.39 mm   | 95.82 mm   | 100.53 mm  |
| Native annulus area (mm <sup>2</sup> )           | 314-380 mm <sup>2</sup>  | 410-500 mm <sup>2</sup>   | 510-630 mm <sup>2</sup>  | 630-770 mm <sup>2</sup>  | 700-840 mm <sup>2</sup>  |
| Area-derived diameter                            | 20-22 mm   | 22.8-25.2 mm  | 25.5-28.3 mm   | 28.3-31.3 mm   | 29.9-32.7 mm   |
| Native annulus size by TEE                       | 17.5-20.5 mm   | 19.5-23.5 mm  | 22.5-26.5 mm   | 25.5-29.5 mm   | 27-31 mm   |

14Fr Python - Expandable Introducer Sheath can be used for all Myval THV diameter sizes up to 32 mm.

\* Contact Meril Representative for information and availability in your country. Myval THV diameter sizes 30.5 mm and 32 mm are currently not CE approved. Myval THV diameter size 30.5 mm is currently pending approval from CDSCO.

# Myval THV: Precise Crimping





## Upon Crimping:

- V-shaped hinges on hexagonal frame fold, generating the dense bands on fluoroscopy
- Vertical "I" connectors give rise to light bands
- Alternating V-folds & vertical connectors give Myval a unique appearance on fluoroscopy for ease of positioning


Open Cells 53%

Closed Cells 47%



# Myval THV: Precise Placement Technique


## Schematic of Myval THV - Ideal Landing Zone



The characteristic bands may not be visible in-case the THV system is not co-axial to the annular plane.  
In this case, THV landing zone must be referenced using Mid-Balloon marker which has to be ~ 3 mm above the annular plane.

## Myval THV: Ground Zero Deployment

- Shallow deployment of Myval with least engagement within LVOT is possible
- Optimal orthotopic anchorage of THV post deployment (80:20) without risk of THV migration
- Minimal infra-annular depth (range of 3-3.50 mm) to avoid conduction system interference (thus minimizing the need of new permanent pacemaker dependency)



Myval THV infra-annular depth range  
3 – 3.50 mm based on THV Ø




# Myval THV Detailed Sizing Guide

|                                 |         |       |       |      |      |       |       |       |      |      |      |       |       |
|---------------------------------|---------|-------|-------|------|------|-------|-------|-------|------|------|------|-------|-------|
| 3D Annular area mm <sup>2</sup> |         | 270   | 280   | 290  | 300  | 310   | 314   | 320   | 330  | 340  | 350  | 363   |       |
| 3D area derived diameter mm     |         | 18.5  | 18.9  | 19.2 | 19.5 | 19.9  | 20.0  | 20.2  | 20.5 | 20.8 | 21.1 | 21.5  |       |
| % Annular area over/under       | 20 mm   | 16.4% | 12.2% | 8.3% | 4.7% | 1.3%  |       |       |      |      |      |       |       |
|                                 | 21.5 mm |       |       |      |      |       | 16%   | 13%   | 10%  | 7%   | 4%   |       |       |
|                                 | 23 mm   |       |       |      |      |       |       |       |      |      |      | 14.5% |       |
| 3D Annular area mm <sup>2</sup> |         | 370   | 380   | 390  | 400  | 410   | 415   | 420   | 430  | 440  | 450  | 460   | 471   |
| 3D area derived diameter mm     |         | 21.7  | 22.0  | 22.3 | 22.6 | 22.8  | 23.0  | 23.1  | 23.4 | 23.7 | 23.9 | 24.2  | 24.5  |
| % Annular area over/under       | 23 mm   | 12.3% | 9.3%  | 6.5% | 3.9% | 1.3%  |       |       |      |      |      |       |       |
|                                 | 24.5 mm |       |       |      |      | 15.0% | 13.6% | 12.2% | 9.6% | 7.1% | 4.8% | 2.5%  |       |
|                                 | 26 mm   |       |       |      |      |       |       |       |      |      |      |       | 12.7% |
| 3D Annular area mm <sup>2</sup> |         | 480   | 490   | 500  | 510  | 520   | 531   | 540   | 550  | 560  | 570  | 580   | 594   |
| 3D area derived diameter mm     |         | 24.7  | 25.0  | 25.2 | 25.5 | 25.7  | 26.0  | 26.2  | 26.5 | 26.7 | 26.9 | 27.2  | 27.5  |
| % Annular area over/under       | 26 mm   | 10.6% | 8.4%  | 6.2% | 4.1% | 2.1%  |       |       |      |      |      |       |       |
|                                 | 27.5 mm |       |       |      |      | 14.2% | 11.9% | 10.0% | 8.0% | 6.1% | 4.2% |       |       |
|                                 | 29 mm   |       |       |      |      |       |       |       |      |      |      | 13.9% | 11.2% |
| 3D Annular area mm <sup>2</sup> |         | 600   | 610   | 620  | 630  | 640   | 650   | 661   | 670  | 680  | 690  | 700   | 710   |
| 3D area derived diameter mm     |         | 27.6  | 27.9  | 28.1 | 28.3 | 28.5  | 28.8  | 29.0  | 29.2 | 29.4 | 29.6 | 29.9  | 30.1  |
| % Annular area over/under       | 29 mm   | 10.1% | 8.3%  | 6.5% | 4.8% |       |       |       |      |      |      |       |       |
|                                 | 30.5 mm |       |       |      |      | 14.2% | 12.4% | 10.5% | 9.0% | 7.4% | 5.9% | 4.4%  |       |
|                                 | 32 mm   |       |       |      |      |       |       |       |      |      |      | 14.9% | 13.3% |
| 3D Annular area mm <sup>2</sup> |         | 720   | 731   | 740  | 750  | 760   | 770   | 780   | 790  | 804  |      |       |       |
| 3D area derived diameter mm     |         | 30.3  | 30.5  | 30.7 | 30.9 | 31.1  | 31.3  | 31.5  | 31.7 | 32.0 |      |       |       |
| % Annular area over/under       | 32 mm   | 11.7% | 10.0% | 8.7% | 7.2% | 5.8%  | 4.4%  | 3.1%  | 1.8% | 0.0% |      |       |       |

Area derived diameter based on MSCT cross-sectional measurements.

# Myval THV Post Deployment Dimension Chart




Largest circumscribable diameter in Open Cell

| Myval Sizes                            | 20 mm    | 23 mm   | 26 mm   | 29 mm    |
|--|----------|---------|---------|----------|
| Myval total height                     | 17.50 mm | 18 mm   | 19 mm   | 20.50 mm |
| Myval open cells 53%                   | 9.20 mm  | 9.46 mm | 9.99 mm | 10.79 mm |
| Myval closed cells 47%                 | 8.15 mm  | 8.39 mm | 8.86 mm | 9.56 mm  |
| Myval Infra-annular height             | 3.05 mm  | 2.85 mm | 3.05 mm | 3.35 mm  |
| Supra-annular height of closed cells   | 5.10 mm  | 5.54 mm | 5.81 mm | 6.21 mm  |
| Recommendation for coronary protection | 10 mm    | 10 mm   | 10 mm   | 10 mm    |

- Consider protection of coronary arteries with a DES especially if height of coronary ostium is < 10 mm from the annular plane and in conjunction with sinus of valsalva dimensions i.e. height & diameters.

# Myval THV Post deployment Dimension Chart: Additional Sizing




Largest circumscribable diameter in Open Cell

| Myval Sizes                            | Myval Intermediate Size |          |          | Myval XL Size |          |
|--|-------------------------|----------|----------|---------------|----------|
|  | 21.5 mm                 | 24.5 mm  | 27.5 mm  | 30.5 mm       | 32 mm    |
| Myval total height                     | 18.50 mm                | 18.90 mm | 19.40 mm | 21.05 mm      | 21.30 mm |
| Myval open cells 53%                   | 9.73 mm                 | 9.94 mm  | 10.20 mm | 11.08 mm      | 11.21 mm |
| Myval closed cells 47%                 | 8.62 mm                 | 8.81 mm  | 9.05 mm  | 9.82 mm       | 9.94 mm  |
| Myval Infra-annular height             | 3.20 mm                 | 2.95 mm  | 3.15 mm  | 3.45 mm       | 3.55 mm  |
| Supra-annular height of closed cells   | 5.42 mm                 | 5.86 mm  | 5.90 mm  | 6.37 mm       | 6.39 mm  |
| Recommendation for coronary protection | 10 mm                   | 10 mm    | 10 mm    | 10 mm         | 10 mm    |

- A balloon occlusion test may be considered to assess the propensity for coronary occlusion. Balloon diameter approximated to shortest axis of CT-derived annular diameter to be considered.

# Navigator THV Delivery System




Navigator – THV Delivery System has been indigenously developed by Meril Life Sciences Pvt. Ltd., India.

# 30.5mm & #32mm dia. are currently not CE approved and approval pending from CDSCO for 30.5mm dia.


% 35mm length is currently not CE approved.

S For Balloon length of 30mm, balloon diameters of 20mm, 21.5mm, 23mm, 24.5mm, 26mm, 27.5mm and 29mm are CE approved.

# Navigator THV Delivery System: Proximal Assembly



Hi-flexion feature ensures tracking the THV delivery system via inner aortic arch curve thereby avoiding contralateral wall scraping.



Caution : Always remember to fully un-flex the Navigator system while withdrawing

# Navigator Balloon Expansion

Navigator balloon with dual expansion ports at each end ensures rapid, simultaneous, controlled expansion (dog-boning) of distal and proximal ends

This typical dog bone pattern of inflation steadies the valve during expansion phase, ensuring its precise annular position and deployment without any risk of valve migration or embolization



Rapid balloon inflation using an inflation device is possible with controlled palm thrust

Rapid balloon deflation within 3-5 sec ensures procedure safety and compliance



# Mammoth - OTW Balloon Catheter

It is not mandatory to pre-dilate the native annulus prior to Myval THV implantation. Operators may consider pre-dilatation based on anatomical considerations.



| Mammoth Balloon Ø  | #14 mm    | 16 mm | 18 mm     | 20 mm | 23 mm | 25 mm | #28 mm | #30 mm |
|--|-----------|-------|-----------|-------|-------|-------|--------|--------|
| Balloon length   | ← 40 mm → |       |           |       |       |       |        |        |
| Balloon RBP  | ← 8 atm → |       | ← 6 atm → |       |       |       |        |        |
| Volume of 75:25 saline:contrast to achieve stated diameter | 8 ml      | 10 ml | 13 ml     | 16 ml | 23 ml | 25 ml | 34 ml  | 42 ml  |

Mammoth – OTW Balloon Catheter has been indigenously developed by Meril Life Sciences Pvt., Ltd. India.

#14 mm, 28 mm, 30 mm dia. are currently not CE approved.  
Balloon RBP is 8 ATM for 14 mm Mammoth Balloon


# Python – Introducer Sheath

## 14Fr Profile for all Myval Ø 20 mm to 32 mm

Sheath expands momentarily like a python swallowing its prey  
 Conveniently allows passage of crimped Myval THV System

14Fr Distal Entry Profile, Allows Atraumatic Percutaneous Access

**Easy retrieval of crimped valve in case of inability to cross annulus**



A separate loading tube ensures temporary opening of hemostatic valves in proximal port allowing smooth passage of crimped Myval THV System

| Common Femoral Artery* Ø (mm) | Myval THV Ø (mm)               |
|-------------------------------|--------------------------------|
| ≥ 5.50 mm                     | 20 mm, 21.5 mm, 23 mm, 24.5 mm |
| ≥ 6.00 mm                     | 26 mm, 27.5 mm, 29 mm          |
| ≥ 6.50 mm                     | 30.5 mm, 32 mm                 |

\*CFA Ø must be MSCT derived. Excluding circumferential Ca<sup>2+</sup>

Python - Introducer Sheath has been indigenously developed by Meril Life Sciences Pvt. Ltd., India.

#Suture eyelets are not yet CE approved.


\$Luer hub design is still not yet CE approved.



# Myval THV Global Clinical Program



2000 Patients



1. Primary endpoint achieved. Presented at PCR LV 2019.

2. 30-days results to be presented at EuroPCR 2020.

3. Initiated. First Patient First Visit expected Apr 2020.

4. Initiated. First Patient First Visit expected Jun 2020.


5. Pre-study activities initiated.

# MyVal-1: Study Design

A prospective, multicentre, single-arm, open-label study of Myval THV in the treatment of severe symptomatic native aortic valve stenosis

Total number of patients: 100

## CLINICAL FOLLOW-UP



|                             |     |     |     |     |
|-----------------------------|-----|-----|-----|-----|
| Clinical Follow-up          | 100 | 100 | 100 | 100 |
| Echocardiographic Follow-up | 100 | 100 | 100 | 100 |
| Quality of Life Measurement | 100 | 100 | 100 | 100 |

Device Sizes – 20, 21.5, 23, 24.5, 26, 27.5 and 29 mm

Dr. Samin Sharma - Chairman, New York - USA

Dr. Ashok Seth - Principal Investigator, New Delhi - India

Dr. Praveen Chandra - Co-ordinating PI, New Delhi - India

Dr. Ravinder Singh Rao - Co-ordinating PI, Jaipur - India

Dr. PK Goel - Scientific Advisor, Lucknow, India

# MyVal-1: Clinical outcomes up to 1-month follow-up

## Excellent clinical safety & efficacy

| Major Adverse Cardiac, Cerebrovascular Renal Events (MACCRE) | Post-procedure | 1-Month Follow-Up |
|--|----------------|-------------------|
| Survival   | 98%            | 97%               |
| Stroke   | 1%             | 2%                |
| Acute renal failure  | 2%             | 2%                |
| Life-threatening or disabling bleeding                       | 1%             | 1%                |
| Myocardial infarction  | 0%             | 0%                |
| Major vascular complications                                 | 1%             | 1%                |
| Minor vascular complications                                 | 2%             | 2%                |
| Repeat hospitalization                                       | NA             | 8%                |
| New permanent pacemaker                                      | 2%*            | 2%                |


\*One patient had RBBB pre-procedure

## MyVal-1: Echocardiographic Findings at 1-month Follow-up

### Echocardiographic findings


| Parameters                                   | Baseline    | Post-procedure | 30-day     |
|--|-------------|----------------|------------|
| Effective orifice area, (cm <sup>2</sup> )   | 0.6 ± 0.2   | 1.7 ± 0.3      | 1.7 ± 0.5  |
| Mean aortic-valve gradient, (mmHg)           | 47.4 ± 8.8  | 8.0 ± 2.7      | 8.8 ± 2.5  |
| Peak aortic-valve gradient, (mmHg)           | 71.7 ± 13.0 | 14.4 ± 2.4     | 15.7 ± 2.8 |
| Trans-aortic velocity, (m/s)                 | 4.5 ± 0.4   | 1.9 ± 0.4      | 1.8 ± 0.4  |
| Mean LVEF, (%)                               | 45.5 ± 11.5 | 47.8 ± 11.1    | 48.6 ± 8.9 |
| Moderate or severe mitral regurgitation, (n) | 2           | 0              | 0          |
| Aortic regurgitation, (n)                    | -           | 0              | 0          |

# Sustained Low Mean Gradients Post-Procedure and ~1.7cm<sup>2</sup> Large EOA at 1-month Follow-up (p<0.0001)



# MyVal-1: Marked improvement in Quality of Life (QoL) parameters


## NYHA Functional Class



| Parameters  | Baseline (n=100) | 1 Month (n=96) |
|---|------------------|----------------|
| Six minute walk test (p<0.0001)                     | 171.7            | 230.69         |
| Kansas city cardiomyopathy Questionnaire (p<0.0001) | 35.27            | 47.91          |

# Interesting Case with Myval THV



71 Years/Male | Symptomatic AS | NYHA class IV | LVEF 25% | Normal Renal Function | Hypertensive, AF on treatment | Moderate Pulmonary Disease | STS Score 7.97% | Valve Type Tri-Leaflet | Peak Velocity 4.59 m/s | Mean gradient 55 mmHg | Peak Gradient 84 mmHg | EOA 0.5 cm<sup>2</sup> | EF 25%



# Myval THV: Sizing Rationale

- Ideal valve sizing may be in the range of 10-15% higher than area derived annulus diameter to have a good valve apposition
- For an area of 561.2mm<sup>2</sup>, a 29mm Myval would be 17.5% over sized
- Considering there is a large Ca<sup>2+</sup> at the LCC running towards LVOT, it was considered to use 2cc less volume to avoid risk of annular rupture

| 3D Annular area mm <sup>2</sup>      |       | 546          | 550          | 560          | 570          | 580          | 590          | 600         | 610         | 615         | 620         | 630         | 640         | 650         | 660          |
|--------------------------------------|-------|--------------|--------------|--------------|--------------|--------------|--------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|--------------|
| 3D area derived diameter mm          |       | 26.4         | 26.5         | 26.7         | 26.9         | 27.2         | 27.4         | 27.6        | 27.9        | 28.0        | 28.1        | 28.3        | 28.5        | 28.8        | 29.0         |
| 29mm<br>% Annular area<br>over/under | 20 mm | -42.7%       | -43.1%       | -44.1%       | -45.1%       | -46.0%       | -46.9%       | -47.8%      | -48.7%      | -49.1%      | -49.5%      | -50.3%      | -51.1%      | -51.8%      | -52.6%       |
|                                      | 23 mm | -24.2%       | -24.7%       | -26.1%       | -27.4%       | -28.6%       | -29.8%       | -31.0%      | -32.1%      | -32.7%      | -33.2%      | -34.3%      | -35.3%      | -36.3%      | -37.3%       |
|                                      | 26 mm | -3.1%        | -3.8%        | -5.5%        | -7.2%        | -8.8%        | -10.3%       | -11.8%      | -13.3%      | -14.0%      | -14.7%      | -16.0%      | -17.3%      | -18.6%      | -19.8%       |
|                                      | 29 mm | <b>20.5%</b> | <b>19.6%</b> | <b>17.5%</b> | <b>15.4%</b> | <b>13.4%</b> | <b>11.5%</b> | <b>9.7%</b> | <b>7.9%</b> | <b>7.0%</b> | <b>6.1%</b> | <b>4.4%</b> | <b>2.8%</b> | <b>1.2%</b> | <b>-0.3%</b> |



Analysis done using CT images and 3mensio software

Clinical case, images and videos courtesy: Dr. Samin Sharma & Dr. Ravinder Singh Rao, EHCC, Jaipur, India




# Case Summary

Baseline - Aortogram



Myval 29mm Positioning



Myval 29mm Final Result



Baseline Echo Results

| Parameters                 | Values              |
|----------------------------|---------------------|
| Valve type                 | Tri-leaflet         |
| Peak Trans-aortic velocity | 4.59 m/s            |
| Mean Trans-aortic gradient | 55 mmHg             |
| Peak gradient              | 84 mmHg             |
| Calculated EOA             | 0.5 cm <sup>2</sup> |
| Severity of AR             | Mild                |
| Severity of MR             | No                  |
| Ejection fraction          | 25%                 |



Post-Procedure Echo Results

| Parameters                 | Values               |
|----------------------------|----------------------|
| Valve type                 | Tri-leaflet          |
| Peak Trans-aortic velocity | 1.20 m/s             |
| Mean Trans-aortic gradient | 3.53 mmHg            |
| Peak gradient              | 5.80 mmHg            |
| Calculated EOA             | >1.0 cm <sup>2</sup> |
| Severity of AR             | No                   |
| Severity of MR             | No                   |
| Ejection fraction          | 25%                  |



# Case Follow-up

| Echo TTE Results at 1 month follow-up | Values               |
|---------------------------------------|----------------------|
| Peak trans-aortic velocity            | 1.76 m/s             |
| Mean trans-aortic gradient            | 6.0 mmHg             |
| Peak gradient                         | 12.0 mmHg            |
| Calculated EOA                        | >1.0 cm <sup>2</sup> |

| Echo TTE Results at 6 months follow-up | Values               |
|--|----------------------|
| Peak trans-aortic velocity             | 1.64 m/s             |
| Mean trans-aortic gradient             | 5.0 mmHg             |
| Peak gradient                          | 11.0 mmHg            |
| Calculated EOA                         | >1.0 cm <sup>2</sup> |



# Myval THV System and Components



### Myval THV Ordering Information

| Diameters    | 20.0 mm | 21.5 mm | 23.0 mm | 24.5 mm | 26.0 mm | 27.5 mm | 29 mm  |
|--------------|---------|---------|---------|---------|---------|---------|--------|
| Product code | MVL200  | MVL215  | MVL230  | MVL245  | MVL260  | MVL275  | MVL290 |

### Navigator - THV Delivery System Ordering Information

| Diameters    | 20.0 x 30 mm | 21.5 x 30 mm | 23.0 x 30 mm | 24.5 x 30 mm | 26.0 x 30 mm | 27.5 x 30 mm | 29 x 30 mm |
|--------------|--------------|--------------|--------------|--------------|--------------|--------------|------------|
| Product code | NVT20030     | NVT21530     | NVT23030     | NVT24530     | NVT26030     | NVT27530     | NVT29030   |

### Mammoth PTA Balloon Ordering Information

| Diameters    | 14 x 40 mm | 16 x 40 mm | 18 x 40 mm | 20 x 40 mm | 23 x 40 mm | 25 x 40 mm | 28 x 40 mm | 30 x 40 mm |
|--------------|------------|------------|------------|------------|------------|------------|------------|------------|
| Product code | MTV1440    | MTV1640    | MTV1840    | MTV2040    | MTV2340    | MTV2540    | MTV2840    | MTV3040    |

### Python 14 Fr Expandable Sheath Ordering Information

| Product code | PHT14 |
|--------------|-------|
|--------------|-------|

### Val-de-Crimp - Heart Valve Crimping Tool

| Product code | VLDC |
|--------------|------|
|--------------|------|

### Atrion QL4015 Inflation Device 15 ATM (40 mL)

| Product code | MVINF15A |
|--------------|----------|
|--------------|----------|

### Syringe Luer Lock 60mL

| Product code | MVSNG60ML |
|--------------|-----------|
|--------------|-----------|

### CrocoDial - Heart Valve Crimping Tool

| Product code | CCD |
|--------------|-----|
|--------------|-----|

### Basix Touch Inflation Device 35 ATM (30 mL)

| Product code | MVINF35A |
|--------------|----------|
|--------------|----------|

### Syringe Luer Lock 50mL

| Product code | MVSNG50ML |
|--------------|-----------|
|--------------|-----------|

#### Manufacturer:

#### Meril Life Sciences Pvt. Ltd.

Survey No. 135/139,  
Bilakhia House, Muktanand Marg,  
Chala, Vapi 396 191.  
Gujarat, India.  
T: +91 260 305 2100  
F: +91 260 305 2125

#### Subsidiary companies:

#### Meril Life Sciences Pvt. Ltd.

301, A-Wing, Business Square,  
Chakala, Andheri Kurla Road,  
Andheri East, Mumbai 400 093  
T +91 22 39350700  
F +91 22 39350777

#### Meril, Inc.

2436 Emrick Boulevard,  
Bethlehem, PA 18020  
T +610 500 2080  
F +610 317 1672

#### Meril South America

Doc Med LTDA  
Al. dos Tupiniquins,  
1079 - Cep: 04077-003 - Moema.  
Sao Paulo, Brazil.  
T +55 11 3624 5935  
F +55 11 3624 5936

#### Meril GmbH.

Bornheimer Strasse 135-137,  
D-53119 Bonn.  
Germany.  
T +49 228 7100 4000  
F +49 228 7100 4001

#### Meril Tibbi Cihazlar

Meril Tibbi Cihazlar İmalat ve Ticaret A.Ş.  
İçerenköy Mah. Çetinkaya Sok.  
Prestij Plaza No:28  
Kat:4 Ataşehir, 34752  
İstanbul / Turkey  
T +90 216 641 44 24  
F +90 216 641 44 25

#### Meril China Co. Ltd.

2301b,23f, Lixin Plaza,  
no 90, South Hubin Road,  
Xiamen, China  
T 0086-592-5368505  
F 0086-592-5368519

#### Meril SA Pty. Ltd.

102, 104, S101 and S102,  
Boulevard West Office Park,  
142 Western Service Road,  
Erf 813 Woodmead Extension 17  
Sandton, Johannesburg – 2191  
South Africa  
T +27 11 465-2049  
F +27 86 471 7941

#### Meril Medical LLC.

Nauchnyi Proezd 19,  
Moscow, Russia – 117 246.  
Office - +7 495 772 7643

#### EU representative.

Obelis S.A.  
Bd, General Wahis 53,  
1030, Brussels, Belgium.  
T +32 2 732 5954  
F +32 2 732 6003  
E mail@obelis.net