

MIRUS™ Skin Stapler

Outstanding performance across wider skin wounds

Skin Staplers are typically used at the end of a surgical procedure for routine wound closure. The staples are then removed with staple extractor after the skin has had an opportunity to heal.

The MIRUS™ Skin Stapler's handle design is characterized by good texture for an user friendly experience.


The angled head provides a clear view to ensure precise staple placement and allows the pins to staple the tissue easily.

Staple material : 316L Stainless Steel
 Sterilization : Ethylene Oxide
 Packaging : 8 pcs/Box
 Code : MSS35 & MSSP35
 Size : 35
 Staples : Regular

Angled head contains 35 staples



Staple Dimension:



Larger diameter wire, greater width, and extra staple height leads to a secure wound closure

Stapler	Skin Stapler
Procedure Application	Wound closure